


## Hématoxyline Harris acidifiée - Haematoxylin Harris Acidified( BMSCR0043-AC-500 )

Libellé	Désignation
Catégorie	COLORANTS/REACTIFS
Code catégorie	CR
Libellé produit	Hématoxyline Harris acidifiée - Haematoxylin Harris Acidified
Code produit	BMSCR0043-AC-500
Image	
Description	<p>Hématoxyline Harris utilise la maturation chimique pour produire une oxydation de l'hématoxyline destinée à une utilisation immédiate dans les laboratoires modernes. Il est souvent utilisé dans les méthodes de coloration régressive, mais peut également être utilisé dans les méthodes progressives en cytologie et comme contre-coloration en histologie.</p> <p style="text-align: center;">L'hématoxyline de Harris colore le noyau d'une couleur bleu-violet foncé.</p> <p style="text-align: center;">Ce produit est la version acidifiée.</p> <p>-----</p> <p>Haematoxylin Harris uses chemical ripening to produce oxidation of haematoxylin for immediate use in modern histology laboratories. It is often used in regressive staining methods, but may also be used progressively in cytology and as a counterstain in histology.</p> <p style="text-align: center;">Harris hematoxylin stains the nuclear with a deep purple-blue colour.</p> <p style="text-align: center;">This is the acidified version.</p>
Prix	22 600 frs