

info@biomedsupexpress.com
 biomedsupexpress@gmail.com
 www.biomedsupexpress.com
 Agrément 4345A/MINSANTE/SG/PDML/SA
 BP: 5084 Nlongkak, Yaoundé - Cameroun (Situé face
 Hôpital Gynéco-Obstétrique de Yaoundé entrée Hôtel
 Segmad)
 (+237) 686 729 251 | 652 505 396



REPUBLIQUE DU CAMEROUN
 Paix - Travail - Patrie

Huile d'Immersion - Immersion Oil(BMSCR0044)

Libellé	Désignation
Catégorie	COLORANTS/REACTIFS
Code catégorie	CR
Libellé produit	Huile d'Immersion - Immersion Oil
Code produit	BMSCR0044
Image	

Libellé	Désignation
Description	Un liquide stable, non séchant et non durcissant, sans danger pour la santé, à utiliser avec des optiques de valeur et des échantillons de lames irremplaçables. Spécification: Apparence: Limpide et exempte de particules Viscosité*: 1 000-1 200 cSt Transmission: 400-900 nm Densité relative*: 1,329 Indice de réfraction*: 1,5182 +/- 0,002 Conserver entre +4 °C et +35 °C *@ 23 °C ----- A stable, non-drying, non-hardening fluid with no health hazard, safe to use with valuable optics and irreplaceable slide specimens. Specification: Appearance: Bright clear and free of particulates Viscosity*: 1000-1200 cSt Transmission: 400-900nm Relative Density*: 1.329 Refractive Index*: 1.5182 +/- 0.002 Store between +4 °C and +35 °C *@ 23 °C
Prix	120 000 frs

Fait à Yaoundé NGOUSSO, le 10-05-2025