

info@biomedsupexpress.com
 biomedsupexpress@gmail.com
 www.biomedsupexpress.com
 Agrément 4345A/MINSANTE/SG/PDML/SA
 BP: 5084 Nlongkak, Yaoundé - Cameroun (Situé face
 Hôpital Gynéco-Obstétrique de Yaoundé entrée Hôtel
 Segmad)
 (+237) 686 729 251 | 652 505 396



REPUBLIQUE DU CAMEROUN
 Paix - Travail - Patrie

Encre de Chine - Indian Ink(BMSCR0045)

Libellé	Désignation
Catégorie	COLORANTS/REACTIFS
Code catégorie	CR
Libellé produit	Encre de Chine - Indian Ink
Code produit	BMSCR0045
Image	

Libellé	Désignation
Description	<p>L'encre de Chine est utilisée en histologie pour marquer de façon permanente les marges des échantillons de tissus, ce qui en fait un outil très utile pour les pathologistes, notamment dans le domaine du diagnostic du cancer.</p> <p>En microbiologie, il est utilisé comme coloration négative, colorant en noir une lame contenant des bactéries afin de voir si une cellule contient une capsule gélatineuse. Son utilisation est recommandée pour l'identification de Cryptococcus neoformans.</p> <p>-----</p> <p>Indian ink is used in histology to permanently mark the margins of tissue specimens, making it a very useful tool to pathologists especially in the cancer diagnosis field.</p> <p>In Microbiology it is used as a negative stain, staining a slide containing bacteria black in order to see if a cell contains a gelatinous capsule. It is recommended for use in the identification of Cryptococcus neoformans.</p>
Prix	120 000 frs

Fait à Yaoundé **NGOUSSO**, le 09-05-2025