

Bromocrésol Violet 0,04% Indicateur alcoolique - Bromocresol Purple 0.04% Alcoholic Indicator(BMSCR0016)

Libellé	Désignation
Catégorie	COLORANTS/REACTIFS
Code catégorie	CR
Libellé produit	Bromocrésol Violet 0,04% Indicateur alcoolique - Bromocresol Purple 0.04% Alcoholic Indicator
Code produit	BMSCR0016
Image	
Description	Bromocrésol violet est un colorant et un indicateur de pH utilisé en laboratoires dans les titrages manuels pour tester la présence de l'albumine dans les échantillons. Sa couleur est jaune quand le pH est en dessous de 5,2 et violet quand il est au-dessus de 6,8. Quand le pH est compris entre 5,2 et 6,8, sa couleur est un mélange de teintes jaunes et violettes. Le changement de couleur est une méthode importante de détection du point final en titrage manuel. ----- Bromocresol violet is a dye and pH indicator used in laboratories in manual titrations to test for the presence of albumin in samples. Its colour is yellow when the pH is below 5.2 and purple when the pH is above 6.8. When the pH is between 5.2 and 6.8, its colour is a gradient mixture of yellow and purple hues. The colour change is the most important method of end point detection in manual titration.
Prix	24 400 frs

Fait à Yaoundé NGOUSSO, le 10-05-2025